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Products for foot- and nail care

# ONLINE SHOP

Instrumente bequem online einkaufen  
**shop.acurata.de**

# INHALT

---

Content

4	<b>VERZÄHNUNGS-NAVIGATOR</b> Toothing-Navigator
6–18	<b>HARTMETALL-INSTRUMENTE</b> Tungsten carbide instruments
6, 18	<b>KUGELFRÄSER</b> Ball cutter
7	<b>HOHLFRÄSER · NAGEL CLEANER</b> Trepine · Nail Cleaner
8	<b>KERAMIKFRÄSER</b> Ceramic cutter
15	<b>HARTMETALLFRÄSER FÜR LINKSHÄNDER</b> Tungsten carbide cutters for sinistrals
16	<b>FRÄSERLEHRE</b> Gage
19–27	<b>DIAMANT-INSTRUMENTE</b> Diamond instruments
28	<b>PODO BASIC SET · NAIL BASIC SET</b>
29	<b>SCHLEIFSTEINE</b> Abrasives
30	<b>SILIKON-POLIERER</b> Silicone-Polishers
31	<b>POLIERBÜRSTEN · ACRYLSTÄNDER</b> Brushes · Acrylic stand
32	<b>INSTRUMENTENSTÄNDER</b> Instrument stands
33	<b>ANWENDUNGS- UND SICHERHEITSHINWEISE</b> General information for Use and Safety
34-35	<b>HYGIENE-EMPFEHLUNG</b> Recommended Hygiene Procedures

---

# VERZÄHNUNGS-NAVIGATOR

Toothing Navigator

## ACURATA HARTMETALLFRÄSER, ROSTFREI

Tungsten carbide cutters, stainless

**NEM** Kombinationsverzahnungen,  
grob und fein (Orange-/Orangering)  
**NEM** Combination cut coarse and fine (orange ring)

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Combination toothing, fine and coarse

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**Characteristic:** fine and coarse segments on one cutter, this means unmatched cutting performance together with smooth surfaces, very smooth running, long service life and especially gentle for the patient  
Applied for German Patent. DE 10 2007 009 304.9

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Quick removal of callus, while smoothing and designing of fungal infected nails; Treatment of cuticle and nail fold.  
**Recommended for diabetics.**

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

Special toothing with cross cut

**Characteristic:**

- All diameters have 8 blades, instead of 6 blades
- Very sharp and effective

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Trephine (Steel instruments, stainless)

	330 104 224 001 018N	330 104 224 001 023N	330 104 224 001 027N
Ø mm / L	1,8 / 6,5	2,3 / 6,5	2,7 / 7,5
€ / VPE	<b>9,90 € / 1</b>	<b>9,90 € / 1</b>	<b>10,90 € / 1</b>

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

- Polygonal working part, without sharp edges
- Gentle removal of debris in the nail fold
- Particularly suitable for extensive treatment

	500 104 189 H 012N	500 104 190 H 012N
Ø mm / L	1,2 / 2,5	1,2 / 5,9
€ / VPE	<b>8,90 € / 1</b>	<b>9,90 € / 1</b>

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**Characteristics:** perfect cutting properties, smooth surfaces, easier cleaning, disinfected and sterilisable, low heat development, smooth running by reduced vibration and especially for risk patients

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- Removal of callus
- Shaping of mycotic nails

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Lamellas toothting with cross cut, very high removal rate

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- Quick removal of thick callus, callosity or thickened nails

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Lamellas toothing, fine with cross cut and PVD-coating (wear protection coating)

**Characteristic:** very high durability, easy cleaning and longer service life, reduced heat development.

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Quick removal of thick grypotic nails and callus, treatment of orthesis silicones and hyperkerathosis, high removal rate

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Quick removal of color gel and ornaments, high removal rate, optimal dust removal

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Lamellas toothing fine with cross cut

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Quick removal of thick gryptic nails and callus, treatment of orthesis silicones and hyperkerathosis, **high removal rate**

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Quick removal of color gel and ornaments, **high removal rate**, optimal dust removal

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Lamellas toothing, coarse with cross cut

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Quick removal of thick gryptic nails and callus treatment of orthesis silicones and hyperkerathosis, **very high removal rate**

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Quick removal of acrylic, color gel and ornaments, very high removal rate

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**Characteristic:** extremely high durability, easy cleaning and longer service life, reduced heat development

	500 104 194 <b>220</b> 031N	500 104 194 <b>220</b> 040N	500 104 137 <b>220</b> 060N	500 104 194 <b>220</b> 050N	500 104 263 <b>220</b> 060N	500 104 274 <b>220</b> 060N
--	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------	-----------------------------

Ø mm / L	3,1 / 13,5	4,0 / 13,5	6,0 / 11,0	5,0 / 14,0	6,0 / 11,5	6,0 / 14,0
----------	------------	------------	------------	------------	------------	------------

€ / VPE	<b>23,90 € / 1</b>	<b>23,90 € / 1</b>	<b>23,90 € / 1</b>	<b>23,90 € / 1</b>	<b>23,90 € / 1</b>	<b>23,90 € / 1</b>
---------	--------------------	--------------------	--------------------	--------------------	--------------------	--------------------

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Delicate treatment of the nail surface, nail wall and cuticle,  
nail pre-polishing

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Removal of callus and thick nails, especially for beginners

	500 104 194 <b>176</b> 040N	500 104 257 <b>176</b> 060N	500 104 274 <b>176</b> 060N	500 104 277 <b>176</b> 040N
Ø mm / L	4,0 / 13,5	6,0 / 12,0	6,0 / 14,0	4,0 / 8,0
€ / VPE	<b>23,90 € / 1</b>	<b>23,90 € / 1</b>	<b>23,90 € / 1</b>	<b>23,90 € / 1</b>

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Careful removal and smoothing of gel, especially for beginners

	500 104 141 <b>134</b> 023N	500 104 194 <b>134</b> 023N	500 104 194 <b>134</b> 040N	500 104 194 <b>134</b> 070N
Ø mm / L	2,3 / 8,0	2,3 / 14,0	4,0 / 12,5	7,0 / 14,0
€ / VPE	<b>14,90 € / 1</b>	<b>18,90 € / 1</b>	<b>23,90 € / 1</b>	<b>23,90 € / 1</b>

	500 104 197 <b>134</b> 014N	500 104 198 <b>134</b> 023N	500 104 274 <b>134</b> 060N	500 104 289 <b>134</b> 023N
Ø mm / L	1,4 / 4,0	2,3 / 8,0	6,0 / 14,0	2,3 / 8,0
€ / VPE	<b>14,90 € / 1</b>	<b>14,90 € / 1</b>	<b>23,90 € / 1</b>	<b>14,90 € / 1</b>

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Treatment of cuticle and nail fold. Smoothing and refinishing

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Non-rounded edges,  
working part completely diamond

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Rounded edges working  
part completely diamond

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Nail smoothening, delicate treatment of cuticle and nails

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Delicate treatment of cuticle prepolishing of nails, final smoothening

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
Polishing of nail prosthesis and nails, medium, grey

	658 104 115 524 060N	658 104 187 524 075N	658 104 200 524 160N	658 104 243 524 055N	658 104 243 524 110N	658 104 273 524 060N	658 104 273 524 100N
Ø mm / L	6,0 / 20,0	7,5 / 20,0	16,0 / 19,0	5,5 / 16,0	11,0 / 20,0	6,0 / 17,0	10,0 / 25,0
€ / VPE	3,50 € / 2	3,50 € / 2	3,50 € / 2	3,50 € / 2	3,50 € / 2	3,50 € / 2	3,50 € / 2

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High-lustre polishing of nail prosthesis and nails, fine, yellow

	658 104 115 514 060N	658 104 187 514 075N	658 104 200 514 160N	658 104 243 514 055N	658 104 243 514 110N	658 104 273 514 060N	658 104 273 514 100N
Ø mm / L	6,0 / 20,0	7,5 / 20,0	16,0 / 19,0	5,5 / 16,0	11,0 / 20,0	6,0 / 17,0	10,0 / 25,0
€ / VPE	3,50 € / 2	3,50 € / 2	3,50 € / 2	3,50 € / 2	3,50 € / 2	3,50 € / 2	3,50 € / 2

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Nylon brush with synthetic filament

**Characteristic:** cleaning is possible by chemical or thermal disinfection

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For storing the cutter at the workplace

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Instrument stand made of aluminium

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
For the safekeeping of rotating instruments at the workplace, during cleaning, disinfection and sterilization. Safe fixing of the instruments by retainer clips.

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Container for 18 handpiece instruments up to an all-over-length of 60 mm

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Instrument stand made of stainless steel

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
Container for 6 handpiece instruments up to an all-over-length of 60 mm

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**37,00 € / 1**

Recommendations for use:

- Use only standard detergents for rotating instruments
- Cleaning in ultrasonic bath is possible
- Sterilizable up to 134 °C
- Do not use the hot air sterilizer and chemiclave because of polymer parts
- Put the instrument stand into the cleaning device so that the instruments are reached by the spray

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General Information for Use and Safety

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acurata instruments are intended for the medical field and must only be used by podiatry or respective experts who, due their training and experience, are familiar with the use of the products. The instruments are supplied non-sterile and have to be prepared prior to first use.

## INFORMATION FOR USE AND SAFETY

Use only technical and hygienic correct drives with a stable ball bearing! Carefully clamp the instruments as deep as possible and assure that the desired speed is reached prior to usage. Please note the maximum allowed engine speed as stated on the packaging label. Avoid levering or canting, otherwise there is risk of fracture. We recommend working with the wet technique; this provides a high level of safety against thermal damages. Separate blunt, twisted or damaged instruments immediately and do not use them again. Avoid contact with H2O2. The contact pressures should be light, do not exceed the maximum of 2N. Excessive contact pressures can cause damage to the instruments. In extreme cases the instrument can fracture! At the same time a higher heat development is generated. Overheating may harm skin / nail, instruments wear faster and surfaces get coarser. Improper usage increases risk, higher wear and inferior work results!

**Remarks:** Polishers and instruments with long working part and respectively long, slender necks or for specific application areas exceeding the maximum allowed speeds can cause oscillations (vibrating) or displacement (whippings). This can result in fracture or heavy damage. Non observance causes an increased security risk! The respective maximum speed is a theoretical value given on the basis of the geometry of the instruments. The choice of the working speed is determined by the material subsequently to be used, the specific indication, the type of handpiece used and the contact pressures and is at users discretion.

### Maximum speeds – table of reference values

**Handpiece instruments (ISO 104)** ISO Ø: 003–023 max. 50.000 min<sup>-1</sup> • 025–040 max. 40.000 min<sup>-1</sup> • 045–060 max. 30.000 min<sup>-1</sup> • 060–250 max. 25.000 min<sup>-1</sup>

**FG-instruments (ISO 314)** ISO Ø: 008–016 max. 450.000 min<sup>-1</sup> • 018–021 max. 300.000 min<sup>-1</sup> • 023–031 max. 160.000 min<sup>-1</sup> • 033–040 max. 120.000 min<sup>-1</sup>

**Polishing instruments** all ISO Ø: max. 5.000–10.000 min<sup>-1</sup>

The information on the label on the packaging applies. Optimum speeds for all instruments depend on the material to be processed. As a general rule 40 – 50 % of the maximum instruction is recommended!

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# RECOMMENDED HYGIENE PROCEDURES

Manufacturer's information on the reprocessing of medical products semi critical B and critical B according to the KRINKO guideline "Requirements on the hygiene for reprocessing of medical products" and the DIN EN ISO 17664.

**Manufacturer:** acurata GmbH & Co. KGaA · Schulstraße 25 · 94169 Thurmansbang · Tel.: +49 8504 9117-0 · Fax.: +49 8504 9117-90 · info@acurata.de · www.acurata.de

**Products:** The present manufacturer's information applies to all instruments supplied by acurata GmbH & Co. KGaA that are used for surgical, periodontal or endodontic treatments. These are rotating tungsten carbide, ceramic, diamond and polishing instruments as well as oscillating instruments made of stainless steel or nickel-titanium. The product range includes exclusively instruments delivered in non-sterile condition. These have to be prepared prior to first use (begin with step 1) and prior to any further use (begin with step 2).

**Limited number of reprocessing cycles:** The end of a product's service life is determined by damage and wear caused by use. Frequent reprocessing does not affect the performance of these instruments. Disposable products (marked  on the label) shall not be reused or reprocessed.

**Basic note:** Observe the legal provisions concerning the reprocessing of medical products valid in your country (e.g. in Germany www.rki.de). The manufacturer assures that the reprocessing methods detailed below are appropriate for the reprocessing of the mentioned groups of instruments for their reutilization in regard of their intended use. The operator is responsible that the applied methods of reprocessing, based on his risk assessment with the used equipment, material, process parameters and staff achieve the required results for the intended use. To guarantee this, routine controls of the validated mechanical and/or manual reprocessing methods are necessary. Any deviation from the validated method below detailed must be checked and released by the operator to ensure effectiveness and to avoid possible adverse consequences. **The hygiene procedure recommendation is available on our website in its current version:**

<https://www.acurata.de/en-US/products/hygiene-recommendations>

## 1 Preparation incl. storage and transportation

For first use preparation begin with step 2. Place instruments immediately or at the latest one hour after use on a patient, in a cleaning/disinfection tank filled with a suitable detergent/disinfectant (non-fixing/ aldehyde-free, e.g. BIB forte eco). Set up the drill bath acc. to manufacturer's instructions; for BIB forte eco mix concentrate with water - first the water, then add concentrate. Cover the tank. Pay attention to the application time (e.g. BIB forte eco 0.5% 60 min.). The transport of the instruments to the place of preparation should be made in a contamination protected cleaning/disinfection tank.

## 2 Cleaning and disinfection

According to the directive of the Robert Koch Institute (RKI) and Commission for Hospital Hygiene and Infection Prevention (KRINKO) it is preferential that the preparation of semi-critical B products is carried out mechanically; critical B products shall in any case be prepared mechanically. For products with long, tight lumina or cavities the cleaning has to be performed mechanical. For root canal instruments silicone stoppers have to be removed prior to reprocessing.

### Mechanical cleaning – validated method

**Equipment:** Cleaning brush, mechanical washer/disinfector acc. to EN ISO 15883 (RDG) (e.g. Miele with Vario TD program), detergent (e.g. 0.5% cleaner Neodisher mediclean), acurata bur stand made of stainless steel.

**Method:** Remove instruments from instrument stand / cleaning-disinfection-tank immediately before mechanical reprocessing and brush off all visible contamination in a cold water bath. Place the instruments into the opened bur stand. Start mechanical cleaning according to the instructions of the manufacturer of the machinery and of the detergent. The following process is validated: Program Vario TD: 2 min. pre-cleaning, 5 min cleaning at 55 °C with detergent, 3 min. neutralizing, 2 min. intermediate rinsing, final rinsing with appropriate (VE-) water 5 min. at >90 °C.

### Thermal disinfection in a validated washer/disinfector

Perform the mechanical cleaning (see above, e.g. Miele washer/disinfector with Vario TD program) incl. thermal disinfection with the instruments fixed in a bur holder. The manufacturer's instructions for the device must be observed. For validated washer/ disinfectors the disinfection is demonstrably assured. acurata products are thermostable up to 134 °C.

### Manual cleaning and disinfection – standardized method

**Equipment:** Cleaning brush (e.g. synthetic brush, sterilizable), ultrasonic bath, detergent & disinfectant with approved efficiency for dental instruments (e.g. BIB forte eco, Alpro Medical), bur stand for rotating oscillating instruments (e.g. acurata bur stand made of stainless steel); the manufacturers' instructions must be observed.

**Method:** Remove instruments from bur stand / cleaning-disinfection tank immediately before manual cleaning. Brush off all visible contamination in a cold water bath. Rinse the instrument and the bur stand under running water. Put the instrument into a suitable strainer element and place it into the ultrasonic unit filled with detergent & disinfectant. Perform cleaning and disinfection according to the instructions of the manufacturers of the ultrasonic bath and the detergent and disinfectant; e.g. BIB forte eco 3% - 10 min. at 55 °C tested to EN 14476. After the application time rinse the instrument thoroughly with appropriate water (e.g. VE-water). Dry the instruments preferably with medical compressed-air. According to KRINKO the manual reprocessing is finalized by a thermal disinfection in a steam sterilizer. Follow manufacturer specifications.

**Visual examination** with a suitable enlarger to ensure that the instrument is clean and undamaged (an enlargement of 8x -10x is recommended). If after reprocessing still residues of contamination are visible, repeat the cleaning and disinfecting process until no visible contamination is left. Instruments showing defects are to be discarded immediately (e.g. missing diamond coating, blunt and chipped blades, deformations, corroded surfaces or non-removable residual contamination).

## 3 Final reprocessing steps – sterile packaging and sterilization – Medical products critical B – validated method with moist heat:


**Equipment:** Steam sterilizer Co. MMM Selectomat HP acurata bur stand made of stainless steel, transparent sterilization bag (Steriking o. VP Stericlin), sealed seam device Co. Hawo

**Packaging:** Prior to sterilization place the instruments in the bur holder and pack them altogether doubly in sterilization bags and weld them with a sealed seam device. The instruments must be protected. For the packing an appropriate standardized method has to be applied.

**Sterilization:** An effective steam sterilization of the packaged instruments is proven successfully in the pre-vacuum steam sterilization method with the following minimal parameters: 3 pre-vacuum phases, 132 °C sterilization temperature, holding time 3 min. (full cycle), drying time 10min. Follow the instructions of the device manufacturer. Note: The products are not suitable for a sterilization in a hot-air sterilizer or chemiclave.

## 4 Transport and storage During transport and storage, the reprocessed products must be protected from recontamination. Further the packed sterile goods must be protected also from dust and moisture.

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acurata GmbH & Co. KGaA · Schulstraße 25 · 94169 Thurmansbang ☎ Telefon +49 8504 9117-0 📠 Fax +49 8504 9117-90