

Crown cutter

FG-grinder for the treatment of ZrO and ceramics



For all dental ceramics, lithium disilicate and zirconium oxide.

A special high performance diamond grain reaches unsurpassed characteristics in cut and durability.



Diamond for ceramics



	L	Ø mm
806 314 139 544 012	5,5	1,2

Diamond for ZrO



	L	Ø mm
806 314 140 526 012	6,5	1,2
806 314 194 526 016	4,5	1,6



Particularities:

- separation and removal of ZrO-crowns
- long lifetime on high performance ceramic
- high cutting performance
- durable
- highly cost efficient
- ideal for turbine, sufficient cooling (50 ml/min⁻¹)


WWW.ACURATA.DE

Certified Quality System DIN EN ISO 9001 · DIN EN ISO 13485

 verkauf@acurata.de

 Schulstraße 25, 94169 Thurmansbang, Germany

 +49 (0) 85 04-91 17 - 15

 +49 (0) 85 04-91 17 - 90