

MicroSIL Lab 95

Blue silicon tray material, addition type, for all laboratory procedures.


Compact material, high precision and simultaneously develops exceptional elasticity



Micropolymer® Srls
dental acrylics production
& distribution

Via Torrette Battifoglia Zona Industriale,
Fraz. S. Andrea delle Fratte,
06132 • Perugia • Italy
Ph. (+39) 0924 49 089
www.micropolymerdental.com
info@micropolymerdental.com

5 kg



Mix Ratio **1 measuring spoon 20g**
• 2,5 cm Cat.

Mixing Time **30"/45"**
Working Time **2'/23°C**
Setting Time **5'**
Final hardness **95sh**

REF 200-800

LOT See on pack

See on pack

DIN Meets the requirements of:
DIN EN 24823 / ISO 4823

CE  Please observe
instruction for use