



HeptaLux

Curing Light

hello



Xpedent Austria GmbH
Millennium Park 9
6890 Lustenau | Austria


andreas@xpedent.com
www.xpedent.com

Distributor :

Pink Hyperspectral LED

- Reduce resin shrinkage, lower the risk of microleakage and secondary caries, and alleviate post-operative sensitivity for patients.
- Enhances the degree of polymerisation, making restorations harder and more wear-resistant.
- Features powerful penetration capabilities; curing depth of up to 8.6 mm.

Normalised intensity (Au) vs Wavelength (nm)




8 modes create clinical flexibility

- Normal Mode
- High Mode
- Super Mode
- Ortho Mode
- Soft Mode
- Pulse Mode
- Check Mode
- Transill Mode
- User memory



Check Mode




Transill Mode

The maximum light intensity reaches 3200mW/cm²



Light intensity measurement



360° rotating head




HeptaLux
11° Illumination angle



Ordinary light guide rod
> 40° Illumination angle

14.5° precise light angle



Clear OLED display Wireless charger

